

# hCG Rapid Pregnancy (Urine) Tests



Vizient Contract LB0812

## Product specifications

- CLIA complexity: Waived for urine
- Cut Off: 20 mIU/mL
- Sensitivities > 99%
- Specificity >99%
- Kit storage condition: Room temperature 39-86°F(2-30°C)
- External controls: Available separately
- Shelf life: 24 months from manufacturer date

## Product Benefits

- Low Cut Off: ability to detect low hCG levels in early pregnancy
- Fast: Provide test result: Urine – 3 minutes
- Sample Type: Urine
- Easy-to-read pink indicator: when sample is well processed
- Reliable: with built-in control ensures accuracy



Item No.	Description	Analytical Sensitivity	Tests/kit	CPT Code
VPT1-CASS	Urine hCG Cassette with Pipette (Novaplust, contract LB0812)	20 mIU/hCG Urine	40	81025
EPT1-CASS	Urine hCG Cassette with Pipette (EKLA)	20 mIU/hCG Urine	40	81025

hCG Rapid Pregnancy Test for the determination of hCG in urine specimens. This test kit is used to obtain a visual qualitative result for early detection of pregnancy.



**EKLA**  
Making a Difference

